## **VIM FISH Probe**

Catalog # FA0247 Size 200 uL

Specification	
Product Description	Made to order FISH probes for identification of gene amplification using Fluorescent In Situ Hybridiz ation Technique. ( <u>Technology</u> ).
Origin	Human
Source	Genomic DNA
Reactivity	Human
Notice	We <b>strongly recommend</b> the customer to use FFPE FISH PreTreatment Kit 1 (Catalog #: <u>KA2375</u> or <u>KA2691</u> ) for the pretreatment of Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections.
Regulation Status	For research use only (RUO)
Supplied Product	DAPI Counterstain (1500 ng/mL ) 250 uL
Storage Instruction	Store at 4°C in the dark.

## Applications

• Fluorescent In Situ Hybridization (Cell)

Protocol Download

Gene Info — VIM	
Entrez GenelD	<u>7431</u>
Gene Name	VIM
Gene Alias	FLJ36605
Gene Description	vimentin

# 😭 Abnova

### **Product Information**

Omim ID	<u>193060</u>
Gene Ontology	Hyperlink
Gene Summary	This gene encodes a member of the intermediate filament family. Intermediate filamentents, along with microtubules and actin microfilaments, make up the cytoskeleton. The protein encoded by thi s gene is responsible for maintaining cell shape, integrity of the cytoplasm, and stabilizing cytoske letal interactions. It is also involved in the immune response, and controls the transport of low-dens ity lipoprotein (LDL)-derived cholesterol from a lysosome to the site of esterification. It functions a s an organizer of a number of critical proteins involved in attachment, migration, and cell signaling. Mutations in this gene causes a dominant, pulverulent cataract
Other Designations	OTTHUMP00000019224

#### Disease

- <u>Alzheimer disease</u>
- Anorexia Nervosa
- <u>Bulimia</u>
- Cognition
- Genetic Predisposition to Disease