

CEN1q (DEAC) FISH Probe

Catalog # FC0009

Size 20 uL

Applications

Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

□

Specification

| | |
|----------------------------|---|
| Product Description | Labeled FISH probes as chromosomal markers using Fluorescent In Situ Hybridization Technique. (Technology). |
| Probe | Name: CEN1q Size: Approximately 600kb Fluorophore: DEAC Location: 1q21.3 |
| Origin | Human |
| Source | Genomic DNA |
| Reactivity | Human |
| Notice | We strongly recommend the customer to use FFPE FISH PreTreatment Kit 1 (Catalog #: KA2375 or KA2691) for the pretreatment of Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections. |
| Regulation Status | For research use only (RUO) |
| Recommend Usage | The probe is provided in 5x concentrated format, to allow mixing of up to 5 chromosome FISH probe s in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer (Cat # U0028 or U0029) before use. |
| Supplied Product | FISH Hybridization Buffer (80 uL) |

Storage Instruction

Store at 4°C in the dark.

Note

Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

Applications

- Fluorescent In Situ Hybridization (Cell)

[Protocol Download](#)