

# CEN9q (DEAC) FISH Probe

Catalog # FC0074

Size 20 uL

## Applications

Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

□

## Specification

|                            |                                                                                                                                                                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Description</b> | Labeled FISH probes as chromosomal markers using Fluorescent In Situ Hybridization Technique. ( <a href="#">Technology</a> ).                                                                                                                                                       |
| <b>Probe</b>               | <b>Name:</b> CEN9q<br><b>Size:</b> Approximately 470kb<br><b>Fluorophore:</b> DEAC<br><b>Location:</b> 9q21                                                                                                                                                                         |
| <b>Origin</b>              | Human                                                                                                                                                                                                                                                                               |
| <b>Source</b>              | Genomic DNA                                                                                                                                                                                                                                                                         |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                                                                               |
| <b>Notice</b>              | We <b>strongly recommend</b> the customer to use FFPE FISH PreTreatment Kit 1 (Catalog #: <a href="#">KA2375</a> or <a href="#">KA2691</a> ) for the pretreatment of Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections.                                                       |
| <b>Regulation Status</b>   | For research use only (RUO)                                                                                                                                                                                                                                                         |
| <b>Recommend Usage</b>     | The probe is provided in 5x concentrated format, to allow mixing of up to 5 chromosome FISH probe s in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer (Cat # <a href="#">U0028</a> or <a href="#">U0029</a> ) before use. |
| <b>Supplied Product</b>    | FISH Hybridization Buffer (80 uL)                                                                                                                                                                                                                                                   |

**Storage Instruction**

Store at 4°C in the dark.

**Note**

Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

## Applications

- Fluorescent In Situ Hybridization (Cell)

[Protocol Download](#)