

CEN12p (Cy5) FISH Probe

Catalog # FC0110 Size 20 uL

Applications

Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

Ш

| Specification | |
|---------------------|---|
| Product Description | Labeled FISH probes as chromosomal markers using Fluorescent In Situ Hybridization Technique. (<u>Technology</u>). |
| Probe | Name: CEN12p Size: Approximately 530kb |
| | Fluorophore: Cy5 Location: 12p11.22 |
| Origin | Human |
| Source | Genomic DNA |
| Reactivity | Human |
| Notice | We strongly recommend the customer to use FFPE FISH PreTreatment Kit 1 (Catalog #: KA2375 or KA2691) for the pretreatment of Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections. |
| Regulation Status | For research use only (RUO) |
| Recommend Usage | The probe is provided in 5x concentrated format, to allow mixing of up to 5 chromosome FISH probe s in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer (Cat # <u>U0028</u> or <u>U0029</u>) before use. |
| Supplied Product | FISH Hybridization Buffer (80 uL) |



Product Information

| Storage Instruction | Store at 4°C in the dark. |
|---------------------|---|
| Note | Hybridization position of the probes on the chromosome: Hybridization position of the probes on the chromosome: |

Applications

• Fluorescent In Situ Hybridization (Cell)

Protocol Download