


Mmp2 (Rat) ELISA Kit

Catalog # KA0393 Size 1 Kit

Applications



The standard curve is for the purpose of illustration only and should not be used to calculate unknowns. A standard curve should be generated each time the assay is performed.

Specification

Product Description	Mmp2 (Rat) ELISA Kit is a sandwich enzyme immunoassay for the quantitative measurement of rat Mmp2.
Suitable Sample	Cell Culture Supernatant, Plasma (Heparin), Serum
Sample Volume	100 uL
Label	HRP-conjugated
Detection Method	Colorimetric
Assay Type	Quantitative
Calibration Range	625 to 40,000 pg/mL
Reactivity	Rat
Regulation Status	For research use only (RUO)
Storage Instruction	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles.

Applications

- Quantification

Gene Info — Mmp2

Entrez GeneID [81686](#)

Gene Name Mmp2

Gene Alias -

Gene Description matrix metallopeptidase 2

Gene Ontology [Hyperlink](#)

Other Designations 72 KDa type IV collagenase|gelatinase A|matrix metalloproteinase 2

Publication Reference

- [Genipin protects against mitochondrial damage of the retinal pigment epithelium under hyperglycemia through the AKT pathway mediated by the miR-4429/JAK2 signaling axis.](#)

Wenshuang Xu, Qingyou Chen, Xiaofeng Zhang, Yue Zhao, Shuang Wu, Chao Yang, Yubao Liu, Lijie Liang, Di Jia, Chaojun Li, Li Fan, Yan Shi.

Annals of Translational Medicine 2022 May; 10(10):587.

- [Intronic regulation of matrix metalloproteinase-2 revealed by in vivo transcriptional analysis in ischemia.](#)

Lee JG, Dahi S, Mahimkar R, Tulloch NL, Alfonso-Jaume MA, Lovett DH, Sarkar R.

PNAS 2005 Nov; 102(45):16345.

- [Structure of human pro-matrix metalloproteinase-2: activation mechanism revealed.](#)

Morgunova E, Tuuttila A, Bergmann U, Isupov M, Lindqvist Y, Schneider G, Tryggvason K.

Science 1999 Jun; 284(5420):1667.

- [Completion of the primary structure of the human type IV collagenase preproenzyme and assignment of the gene \(CLG4\) to the q21 region of chromosome 16.](#)

Huhtala P, Eddy RL, Fan YS, Byers MG, Shows TB, Tryggvason K.

Genomics 1990 Mar; 6(3):554.