

Global Di-Methyl Histone H4K20 Quantification Kit (Colorimetric)

Catalog # KA1519

Size 1 Kit

Specification

Product Description	Global Di-Methyl Histone H4K20 Quantification Kit (Colorimetric) is used for measuring global di-methylation of histone H4K20.
Suitable Sample	Histone Extract
Sample Volume	50 uL (1-2 ug included)
Detection Method	Colorimetric
Calibration Range	20 to 2000 ng/well
Limit of Detection	5 ng/well
Regulation Status	For research use only (RUO)
Storage Instruction	Store the standard at -20°C. Store all other components at 4°C away from light.

Applications

- Quantification