Global Acetyl Histone H4K12 Quantification Kit (Colorimetric)

Catalog # KA1536 Size 1 Kit

Specification	
Product Description	Global Acetyl Histone H4K12 Quantification Kit (Colorimetric) is used for measuring global acetylatio n of histone H4K12.
Suitable Sample	Histone Extract
Sample Volume	50 uL (1-2 ug included)
Detection Method	Colorimetric
Calibration Range	20 to 5000 ng/well
Limit of Detection	2 ng/well
Specificity	Histone extracts from mammalian cells/ tissues
Regulation Status	For research use only (RUO)
Storage Instruction	Store the detection antibody and the standard at -20°C. Store all other components at 4°C away from light.

Applications

Quantification