

Total Histone H4 Acetylation Detection Fast Kit (Colorimetric)

Catalog # KA1540

Size 1 Kit

Specification

Product Description	Total Histone H4 Acetylation Detection Fast Kit (Colorimetric) is used for measuring total acetylation of histone H4.
Suitable Sample	Histone Extract
Sample Volume	50-200 ng
Detection Method	Colorimetric
Calibration Range	5 to 2000 ng/well
Limit of Detection	2 ng/well
Specificity	Histone extracts from mammalian cells/ tissues
Regulation Status	For research use only (RUO)
Storage Instruction	Store the Detection Antibody and the Standard Control at -20°C. Store Stop Solution at room temperature away from light. Store all other components at 4°C away from light.

Applications

- Functional Study