


# Protease Activity Assay Kit

Catalog # KA6048

Size 1 Kit

## Applications



### Sample Data:

The data shown here is an example of fluorescence polarization (FP) or total fluorescence data typically produced with this kit; however, your results will not be identical to these. Do not use the data below to directly compare to your samples. Your results could differ substantially.

## Specification

<b>Product Description</b>	Protease Activity Assay Kit is used for the measurement of protease activity.
<b>Sample Volume</b>	100 uL
<b>Detection Method</b>	Fluorometric
<b>Excitation</b>	480-495 nm
<b>Emission</b>	515-525 nm
<b>Regulatory Status</b>	For research use only (RUO)

**Quality Control Testing**

Digestion of FITC-casein by trypsin, chymotrypsin and proteinase K, (A) Fluorescence polarization, (B) Total fluorescence.

Sample Data:

The data shown here is an example of fluorescence polarization (FP) or total fluorescence data typically produced with this kit; however, your results will not be identical to these. Do not use the data below to directly compare to your samples. Your results could differ substantially.

**Storage Instruction**

Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

## Applications

- Functional Study