


5-Methylcytosine monoclonal antibody, clone 5MC-CD

Catalog # MAB6766 Size 100 ug

Applications




Immunofluorescence


Methylation of chloroplast DNA visualized by immunochemistry. Samples are Chlamydomonas me-1 cells.

Left: DAPI-stained cells.

Middle: Cells stained with 5-Methylcytosine monoclonal antibody, clone 5MC-CD (Cat # MAB6766) and FITC-conjugated 2nd antibody.

Right: Merged image Chloroplast DNA is exclusively methylated in gamete cells.





Immunofluorescence

Detection of DNA methylation in mouse embryonic stem cells by immunofluorescence staining with the 5-Methylcytosine monoclonal antibody, clone 5MC-CD (Cat # MAB6766). Intense 5-methylcytosine staining at pericentromeric regions was seen in the mitotic chromosome and interphase nuclei of ESCs.

| | | | • |
|---|-------|------|-----|
| | pecif | ıcat | ınn |
| U | | IGAL | |
| | | | |

Product Description Mouse monoclonal antibody raised against 5-Methylcytosine.

Immunogen 5-Methylcytosine conjugated with BSA.



Product Information

| Host | Mouse | |
|---------------------|--|--|
| Form | Liquid | |
| Isotype | lgM | |
| Recommend Usage | Immunocytochemistry (1:50-1:100) | |
| | Western Blot (~1:1000) | |
| | The optimal working dilution should be determined by the end user. | |
| Storage Buffer | In PBS (50% glycerol). | |
| Storage Instruction | Store at -20°C. | |
| | Aliquot to avoid repeated freezing and thawing. | |
| | | |

Applications

- Immunocytochemistry
- Immunofluorescence

Methylation of chloroplast DNA visualized by immunochemistry. Samples are Chlamydomonas me-1 cells.

Left: DAPI-stained cells.

Middle: Cells stained with 5-Methylcytosine monoclonal antibody, clone 5MC-CD (Cat # MAB6766) and FITC-conjugated 2nd antibody.

Right: Merged image Chloroplast DNA is exclusively methylated in gamete cells.

Immunofluorescence

Detection of DNA methylation in mouse embryonic stem cells by immunofluorescence staining with the 5-Methylcytosine monoclonal antibody, clone 5MC-CD (Cat # MAB6766). Intense 5-methylcytosine staining at pericentromeric regions was seen in the mitotic chromosome and interphase nuclei of ESCs.

Publication Reference

The SRA protein Np95 mediates epigenetic inheritance by recruiting Dnmt1 to methylated DNA.

Sharif J, Muto M, Takebayashi S, Suetake I, Iwamatsu A, Endo TA, Shinga J, Mizutani-Koseki Y, Toyoda T, Okamura K, Tajima S, Mitsuya K, Okano M, Koseki H.

Nature 2007 Dec; 450(7171):908.

Application: IF, Human, Human embryonic stem cells





• A chloroplast-resident DNA methyltransferase is responsible for hypermethylation of chloroplast genes in Chlamydomonas maternal gametes.

Nishiyama R, Ito M, Yamaguchi Y, Koizumi N, Sano H.

PNAS 2002 Apr; 99(9):5925.

Application: IF, Chlamydomonas reinhardtii, me-1 cells

 Detection of heavy methylation in human repetitive DNA subsets by a monoclonal antibody against 5methylcytosine.

Sano H, Imokawa M, Sager R.

Biochimica et Biophysica Acta 1988 Nov; 951(1):157.

Application: AFC, IA, Bovine, Human, Bovine nuclear DNA, Human lymphocytes DNA