

## Bone Marrow (diabetes) FFPE tissue slide

Catalog # S1174 Size 5 Pieces/Box

Specification	
Product Description	Bone Marrow (diabetes) FFPE tissue slide contains tissue sections from a single human donor.
Tissue Type	Disease
Clinical Diagnosis	Diabetic diseased
Tissue	Bone marrow
Host	Human
Form	Tissue sections are 5 um in thickness, and mounted on positively charged glass slides.
Quality Control Testing	H&E staining was performed for at least one of the tissue slides from each lot.
Recommend Usage	Please bake slides at 60°C for 30 minutes before use.
Storage Instruction	Store at 4°C.
Note	The product is based on the sample repository network established following the IRB-approved ethic al standard and procedures.
Recommend Reagents	IHC antibodies
	<u>CISH probes</u>
	FISH probes
	<u>mutaFISH™ probes</u>
Videos	H&E Staining
	<u>Immunohistochemistry</u>



## **Product Information**

Services FFPE Tumor Tissue DNA Service

Immunohistochemistry Service

Custom FISH Probe Service

Custom mutaFISH™ Probe Service

## **Applications**

- Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)
- Protein Profiling