

NTP Set (100 mM)

Catalog # U0438

Size 1 Set

Specification

Product Description	The NTP Set consists of 4 separate 100 mM solutions of ATP, UTP, CTP, and GTP (pH 7.4) as sodium salts.
Form	Liquid
Purity	>= 98% purity of each component confirmed by HPLC.
Storage Buffer	In 100 mM aqueous solution, pH7.4.
Storage Instruction	Store at -20°C. Avoid repeated freezing and thawing.
Note	Mix thoroughly before use and dilute with nuclease-free water if necessary.

Applications

- In vitro Transcription