


2-NBDG Glucose Uptake Assay Kit

Catalog # KA6077

Size 1 Kit

Applications



Flow Cytometry

Flow cytometry of 2-NBDG uptake in CHO-K1 cells using the 2-NBDG Glucose Uptake Assay Kit. CHO-K1 cells were treated with or without 100 μ M Phloretin at 37°C for 1 hour, then incubated with 100 μ M 2-NBDG staining solution for 20 minutes. To prepare adherent CHO-K1 cells for flow cytometry, EDTA was used to detach cells after staining. Fluorescence intensity was measured using ACEA NovoCyte flow cytometer in FITC channel.

Example Data Analysis and Figures.

Example Data Analysis and Figures.

Fluorescence images of 2-NBDG uptake in CHO-K1 cells using the 2-NBDG Glucose Uptake Assay Kit. CHO-K1 cells at 40,000 cells/well/100 μ L were seeded overnight in a 96-well black wall/clear bottom plate. Cells were treated with 20 mM Glucose (B) or 100 μ M Phloretin (C) at 37°C for 1 hour, then incubated with 100 μ M 2-NBDG staining solution for 20 minutes. Untreated control cells were stained under the same conditions. The fluorescence signal was measured using a fluorescence microscope with FITC filter.

Specification

Product Description

2-NBDG Glucose Uptake Assay Kit provides a sensitive and non-radioactive assay for measuring glucose uptake in cultured cells. The fluorescence signal can be monitored by fluorescence microscope or flow cytometer with a 488 nm laser and 530/30 nm emission filter (FITC channel). This Assay Kit is the most robust tool for monitoring glucose transporters.

Suitable Sample	Cultured Cells.
Detection Method	Fluorimetric
Platform	Instrument: Flow cytometer Excitation: 488 nm laser Emission: 530/30 nm filter Instrument specification(s): FITC channel Instrument: Fluorescence microscope Excitation: FITC filter Emission: FITC filter Recommended plate: Black wall/clear bottom
Regulatory Status	For research use only (RUO)
Storage Instruction	Store the kit at -20°C and avoid from light.
Note	Example Data Analysis and Figures. Example Data Analysis and Figures. Fluorescence images of 2-NBDG uptake in CHO-K1 cells using the 2-NBDG Glucose Uptake Assay Kit. CHO-K1 cells at 40,000 cells/well/100 uL were seeded overnight in a 96-well black wall/clear bottom plate. Cells were treated with 20 mM Glucose (B) or 100 uM Phloretin (C) at 37°C for 1 hour, then incubated with 100 uM 2-NBDG staining solution for 20 minutes. Untreated control cells were stained under the same conditions. The fluorescence signal was measured using a fluorescence microscope with FITC filter.

Applications

- Functional Study
- Flow Cytometry

Flow cytometry of 2-NBDG uptake in CHO-K1 cells using the 2-NBDG Glucose Uptake Assay Kit. CHO-K1 cells were treated with or without 100 uM Phloretin at 37°C for 1 hour, then incubated with 100 uM 2-NBDG staining solution for 20 minutes. To prepare adherent CHO-K1 cells for flow cytometry, EDTA was used to detach cells after staining. Fluorescence intensity was measured using ACEA NovoCyte flow cytometer in FITC channel.

Publication Reference

- [Respiratory syncytial virus co-opts hypoxia-inducible factor-1 \$\alpha\$ -mediated glycolysis to favor the production of infectious virus.](#)

Li-Feng Chen, Jun-Xing Cai, Jing-Jing Zhang, Yu-Jun Tang, Jia-Yi Chen, Si Xiong, Yao-Lan Li, Hong Zhang, Zhong Liu, Man-Mei Li.

mBio 2023 Oct; 14(5):e0211023.

Application: Func, Human, Hep-2 cells

- [SNHG15 promotes chemoresistance and glycolysis in colorectal cancer.](#)

Min Li, Shengbai Sun, Zehua Bian, Surui Yao, Meng Liu, Xiaohong You, Min Li.

Pathology, Research and Practice 2023 Jun; 246:154480.

Application: Func, Human, DLD1, HTC8 cells