

# Protein Staining Solution

Catalog # TM0008

Size 500 mL

## Specification

**Product Description**

Protein Staining Solution stains proteins in polyacrylamide gel after SDS gel electrophoresis (detection limit  $\geq 8$  ng of protein) with \*CoomassieR-G250. All processes including washing and destaining processes can be performed in 1.5 hr. After staining of gel (20 min), proteins can be seen before destaining. Overnight staining and overnight destaining is allowable.

**Storage Instruction**

Store at 4°C.

**Note**

Protein Staining Solution can't be re-used. Since Protein Staining Solution dyes the skin, please wear gloves throughout using Protein Staining Solution.

## Applications

- Electrophoresis