

CEN9q (DEAC) FISH Probe

Catalog # FC0074

Size 20 uL

Applications

Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

□

Specification

Product Description Labeled FISH probes as chromosomal markers using Fluorescent In Situ Hybridization Technique. ([Technology](#)).

Probe
Name: CEN9q
Size: Approximately 470kb
Fluorophore: DEAC
Location: 9q21

Origin Human

Source Genomic DNA

Reactivity Human

Notice We **strongly recommend** the customer to use FFPE FISH PreTreatment Kit 1 (Catalog #: [KA2375](#) or [KA2691](#)) for the pretreatment of Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections.

Regulation Status For research use only (RUO)

Recommend Usage The probe is provided in 5x concentrated format, to allow mixing of up to 5 chromosome FISH probe s in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer (Cat # [U0028](#) or [U0029](#)) before use.

Supplied Product FISH Hybridization Buffer (80 uL)

Storage Instruction

Store at 4°C in the dark.

Note

Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

Applications

- Fluorescent In Situ Hybridization (Cell)

[Protocol Download](#)