

# CEN11q (FITC) FISH Probe

Catalog # FC0101 Size 20 uL

### **Applications**

#### Hybridization position of the probes on the chromosome:

Hybridization position of the probes on the chromosome:

Ш

Specification	
Product Description	Labeled FISH probes as chromosomal markers using Fluorescent In Situ Hybridization Technique. ( <u>Technology</u> ).
Probe	Name: CEN11q Size: Approximately 680kb Fluorophore: FITC Location: 11q12.1
Origin	Human
Source	Genomic DNA
Reactivity	Human
Notice	We <b>strongly recommend</b> the customer to use FFPE FISH PreTreatment Kit 1 (Catalog #: KA2375 or KA2691) for the pretreatment of Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections.
Regulation Status	For research use only (RUO)
Recommend Usage	The probe is provided in 5x concentrated format, to allow mixing of up to 5 chromosome FISH probe s in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer (Cat # <u>U0028</u> or <u>U0029</u> ) before use.
Supplied Product	FISH Hybridization Buffer (80 uL)



### **Product Information**

Storage Instruction	Store at 4°C in the dark.
Note	Hybridization position of the probes on the chromosome: Hybridization position of the probes on the chromosome:

## Applications

• Fluorescent In Situ Hybridization (Cell)

Protocol Download