AMPK Signaling Phospho-Specific Array

Catalog # AA0014 Size 1 Kit (2 Samples)

| Specification | |
|---------------------|--|
| Product Description | AMPK Signaling Phospho-Specific Array includes 174 highly specific and well-characterized phosph orylation antibodies in the AMPK signaling pathway. |
| Substrate | The arrays are printed on standard-size coated glass microscope slides and can be scanned on all microarray scanners that are compatible with 76 x 25 x 1 mm (3 in. x 1 in. x 1mm) slides. |
| Suitable Sample | Cultured Cells, Frozen Tissue, FFPE Tissue, Serum |
| Detection Method | Fluorometric |
| GAL file | GenePix Array List File Download |
| Storage Instruction | Store at 4°C. |
| Note | This product must be accompanied with Antibody Array Detection Kit (<u>AA0135</u>) for use. |

Applications

• Protein Profiling