

## GPCR Signaling to MAPK/ERk Phospho-Specific Array

Catalog # AA0027 Size 1 Kit (2 Samples)

Specification	
Product Description	GPCR Signaling to MAPK/ERk Phospho-Specific Array includes 193 highly specific and well-charac terized phosphorylation antibodies related to this signaling pathway.
Substrate	The arrays are printed on standard-size coated glass microscope slides and can be scanned on all microarray scanners that are compatible with 76 x 25 x 1 mm (3 in. x 1 in. x 1mm) slides.
Suitable Sample	Cultured Cells, Frozen Tissue, FFPE Tissue, Serum
Detection Method	Fluorometric
GAL file	GenePix Array List File  Download
Storage Instruction	Store at 4°C.
Note	This product must be accompanied with Antibody Array Detection Kit (AA0135) for use.

## Applications

Protein Profiling