

100 Base-Pair DNA Ladder, Express

Catalog # D0028

Size 500 uL

Specification

Product Description

Consists of multiple repeats of 100-bp fragments from 100 bp up to 1000 bp. Supplied in sufficient quantity for 100 loadings using 4-5 ul per loading. Not designed for quantitating DNA concentration in a sample. Ladders are available premixed in as ready-to-use format premixed with loading-dye-solution (Express).

Recommend Usage

4 to 5 uL/lane

Storage Instruction

Store at -20°C. Stable at room temperature for 6 months.

Applications

- Electrophoresis