


# Glucose-6-Phosphate Assay Kit (Colorimetric)

Catalog # KA4570

Size 1 Kit

## Applications



The standard curve is for the purpose of demonstration only and should not be used to calculate unknowns. A standard curve should be generated each time the assay is performed.

## Specification

### Product Description

Glucose-6-Phosphate Assay Kit (Colorimetric) is used for quantitative determination of Glucose-6-P phosphate (G6P) using colourmetric methods.

### Suitable Sample

Biological samples (e.g. Plasma, Serum, Tissue and Culture Media.).

### Sample Volume

20 uL

### Detection Method

Colorimetric

### Assay Type

Quantitative

### Calibration Range

10 to 1000 uM

### Regulation Status

For research use only (RUO)

### Quality Control Testing

Standard Curve

The standard curve is for the purpose of demonstration only and should not be used to calculate unknowns. A standard curve should be generated each time the assay is performed.

### Storage Instruction

Store the kit at -20°C.

## Applications

- Functional Study