## miRNA Isolation Kit

Catalog # KA6533 Size 1 Kit

| Specification       |  |
|---------------------|--|
| Product Description | miRNA lsolation Kit is a miRNA extraction and isolation kit based on the adsorption column method. I<br>t can efficiently extract small RNAs such as miRNA, siRNA, snRNA, etc., which are 20-200 nt in lengt<br>h, and can also extract total RNA. |
| Supplied Product    | Kit content:   |
|                     | <ol> <li>Buffer MZ</li> <li>Buffer RW</li> <li>Buffer MRD</li> <li>Rnase-Free ddH<sub>2</sub>O</li> <li>Rnase-Free Columns miRspin set</li> </ol>  |
|                     | 6. Rnase-Free Columns miRelute set<br>7. Rnase-Free Tubes  |
| Pack Size           | 50 reactions   |
| Suitable Sample     | cells, animal tissues, plant tissues, serum, plasma  |
| Sample Size         | 30 - 50 mg of animal tissue, 100 mg plant tissue, 1 x 10 <sup>7</sup> cultured cells   |
| Regulation Status   | For research use only (RUO)  |
| Storage Instruction | Store the kit at room temperature.<br>For best results, please store the kit components as described in the protocol.  |

## Applications

• Purification