

Bioactive

Collagenase Type 2 (C. histolyticum)

Catalog # P5277 Size 50 mg

Specification	
Product Description	Lyophilized. Dialyzed. 0.22um filtered. Contains higher clostripain activity.
Biological function	Collagenases are enzymes that break the peptide bonds in collagen.
Host	Clostridium Species
Form	Lyophilized
Preparation Method	Native protein purified from Clostridium histolyticum
Activity	>=125 units per mg dry weight. One Unit releases one micromole of L-leucine equivalents from colla gen in 5 hours at 37°C, pH 7.5.
Recommend Usage	Collagenase Assay; Bone, heart, liver, thyroid and salivary primary cell isolation.
Storage Instruction	Store at 4°C on dry atmosphere. For primary cell isolation, culture and tissue dissociation: Reconstituted in balanced salt solutions (such as Hank's, Earle's and others) at concentrations from 0.05% to 0.5% (w/v). For enzymetic assay: Reconstituted in 0.05 M TES buffer with 0.36 mM calcium chloride, pH 7.5 at a concentration of 1.0 mg/ml.

Applications

Enzyme Activity