

6 X DNA loading Dye

Catalog # D0008

Size 4 mL

Specification

Product Description	6 X DNA loading Dye is useful for loading DNA samples into wells on agarose gel electrophoresis. It contains EDTA, glycerol, bromphenol blue and 10 mM Tris-HCl (pH8.0).
Recommend Usage	Mix 6 X BPB loading Dye with 1/5 volume of DNA solution to a well on agarose gel electrophoresis.
Storage Instruction	Store at room temperature. For long term storage, store at -20°C.

Applications

- Polymerase Chain Reaction
- Electrophoresis