

# MAPK Phospho-Specific Array

Catalog # AA0033

Size 1 Kit (2 Samples)

## Specification

<b>Product Description</b>	MAPK Phospho-Specific Array includes 185 highly specific and well-characterized phosphorylation antibodies in the MAPK pathway.
<b>Substrate</b>	The arrays are printed on standard-size coated glass microscope slides and can be scanned on all microarray scanners that are compatible with 76 x 25 x 1 mm (3 in. x 1 in. x 1mm) slides.
<b>Suitable Sample</b>	Cultured Cells, Frozen Tissue, FFPE Tissue, Serum
<b>Detection Method</b>	Fluorometric
<b>GAL file</b>	GenePix Array List File <a href="#">Download</a>
<b>Storage Instruction</b>	Store at 4°C.
<b>Note</b>	This product must be accompanied with Antibody Array Detection Kit ( <a href="#">AA0135</a> ) for use.

## Applications

- Protein Profiling