

# FFPE FISH PreTreatment Kit 2

Catalog # KA2376

Size 1 Kit

## Specification

### Product Description

Pre-treatment is an essential step for effective FISH on Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections. This kit contains key reagents used to pretreat FFPE tissue sections from dewaxing with xylene through proteolytic digestion to help permeabilization of cell membranes to facilitate penetration of fluorescence labeled probes. After proteolytic digestion, proceed with subsequent steps with the FISH assay protocol.

### Regulation Status

For research use only (RUO)

### Supplied Product

Kit content:

1. **10x Paraffin Pretreatment Solution**, 500 mL
2. **Protease (Pepsin)**, 0.5 mL x 10 vials

### Storage Instruction

After arrival, please store the **10x Paraffin Pretreatment Solution** at room temperature (RT), and store the **Protease (Pepsin)** at -20°C.

## Applications

- Fluorescent In Situ Hybridization (Cell)

[Protocol Download](#)