


# sgRNA Screening Kit

Catalog # KA6447      Size 1 Kit

## Applications

### Electrophoresis



Electrophoresis efficiency screening of gRNA efficiency screening. (1) gRNA9 without Cas9 (2) Cas9 without gRNA9 (3) gRNA1+Cas9 (4) gRNA2 + Cas9 (5) gRNA3 + Cas9 (6) gRNA4 + Cas9 (7) gRNA5 + Cas9 (8) gRNA6 + Cas9 (9) gRNA7 + Cas9 (10) gRNA8 + Cas9 (11) gRNA9 + Cas9 (12) gRNA10 + Cas9 (13) gRNA 11 + Cas9.

## Specification

<b>Product Description</b>	Assessing sgRNA efficiency before cell transduction, allows you to identify the highly effective CRISPR sgRNA.
<b>Recommend Usage</b>	CRISPR Genome editing The optimal working dilution should be determined by the end user.
<b>Storage Instruction</b>	Store at -80°C. Aliquot to avoid repeated freezing and thawing.

## Applications

- Electrophoresis

Electrophoresis efficiency screening of gRNA efficiency screening. (1) gRNA9 without Cas9 (2) Cas9 without gRNA9 (3) gRNA1+Cas9 (4) gRNA2 + Cas9 (5) gRNA3 + Cas9 (6) gRNA4 + Cas9 (7) gRNA5 + Cas9 (8) gRNA6 + Cas9 (9) gRNA7 + Cas9 (10) gRNA8 + Cas9 (11) gRNA9 + Cas9 (12) gRNA10 + Cas9 (13) gRNA 11 + Cas9.

- CRISPR Genomic Editing