

Proteoliposome Protein Mock (Negative) Control

Catalog # P3689

Size 2 x 35 uL

Specification

Product Description	This is the mock (negative) control of proteoliposome protein prepared using empty expression vect or (without insert).
Form	Liquid
Preparation Method	<i>in vitro</i> wheat germ expression system with proprietary liposome technology
Recommend Usage	We recommend you start with 10 uL Proteoliposome Protein Negative (Mock) Control and 10 uL targ et protein. But the end user should determine the optimal condition.
Storage Buffer	25 mM Tris-HCl of pH8.0 containing 2% glycerol.
Storage Instruction	Store at -80°C. Aliquot to avoid repeated freezing and thawing.
Note	Best use within three months from the date of receipt of this product.

Applications

- Control