


# 3D Magnet Separator


Catalog # DU0009 Size 1 Set

## **Applications**



#### **WBC** Depletion

Representative images of white blood cells before and after WBC depletion using 3D Magnet Separator & CD45 Depletion Kit (Cat# DM0006).



Up to 2.7 log depletion of Jurkat cells

## **Specification**

**Product Description** 3D Magnet Separator contains 88 magnents for isolating magnetic beads.

Amenable Platform LiquidCell™ negative enrichment cell isolation



#### **Product Information**

**Specification** Height: 46 mm

Weight: 86 g

Outer diameter: 51 mm

Inner diameter: 28.4 mm (Max.)

Well: v-shaped bottom

**Product Image** 

Quality Control Testing

Cell Depletion Assay (IF)

Up to 2.7 log depletion of Jurkat cells

Regulatory Status GMP grade

Storage Instruction Store at room temperature

## **Applications**

WBC Depletion

Representative images of white blood cells before and after WBC depletion using 3D Magnet Separator & CD45 Depletion Kit (Cat# DM0006).