

NanoAb™

RecomAb™

## FluoTag-X4 anti-GFP, clone 1H1/1B2 (AF568)

Catalog # RAB00732 Size 200 uL

Specification	
Product Description	FluoTag-X4 anti-GFP two single-domain antibody (sdAb, also known as NanoAb <sup>TM</sup> ) raised against t wo distinct epitopes on GFP and its most common derivatives.
Antibody Species	Camelid
Produced in	E. coli
Specificity	Recognizes GFP (green fluorescent protein) and common GFP derivatives like EGFP, mEGFP, Siri us, tSapphire, Cerulean, eCFP, mTurquoise, acGFP, Emerald, superecliptic pHluorin, paGFP, super folder GFP, eYFP, mVenus and Citrine.
Form	Lyophilized
Conjugation	AZDye 568 Two fluorophores of FluoTag-X4 series are site-specifically coupled to each individual FluoTag mole cule.
Concentration	5 uM
Recommend Usage	Immunofluorescence (1:500) The optimal working dilution should be determined by the end user.
Storage Buffer	Lyophilized from PBS pH 7.4 containing 2.5 uM fluorescently labeled clone 1H1 and 2.5 uM fluoresce ntly labeled clone 1B2. Reconsitute with 200 uL 50% glycerol.
Storage Instruction	Store at -20°C. For long term storage store at -80°C. Avoid prolonged exposure to light. Aliquot to avoid repeated freezing and thawing.
Note	It is selected from the alpaca in vitro synthetic gene library, produced by phage display derived reco mbinant antibodies.

## Applications

Immunofluorescence