

NanoAb™

RecomAb™

FluoTag-X4 anti-GFP, clone 1H1/1B2 (Abberior STAR 635P)

Catalog # RAB00734 Size 200 uL

Specification	
Product Description	FluoTag-X4 anti-GFP two single-domain antibody (sdAb, also known as NanoAb TM) raised against t wo distinct epitopes on GFP and its most common derivatives.
Antibody Species	Camelid
Produced in	E. coli
Specificity	Recognizes GFP (green fluorescent protein) and common GFP derivatives like EGFP, mEGFP, Siri us, tSapphire, Cerulean, eCFP, mTurquoise, acGFP, Emerald, superecliptic pHluorin, paGFP, super folder GFP, eYFP, mVenus and Citrine.
Form	Lyophilized
Conjugation	Abberior STAR 635P Two fluorophores of FluoTag-X4 series are site-specifically coupled to each individual FluoTag mole cule.
Concentration	5 uM
Recommend Usage	Immunofluorescence (1:500) The optimal working dilution should be determined by the end user.
Storage Buffer	Lyophilized from PBS pH 7.4 containing 2.5 uM fluorescently labeled clone 1H1 and 2.5 uM fluorescently labeled clone 1B2. Reconsitute with 200 uL 50% glycerol.
Storage Instruction	Store at -20°C. For long term storage store at -80°C. Avoid prolonged exposure to light. Aliquot to avoid repeated freezing and thawing.
Note	It is selected from the alpaca in vitro synthetic gene library, produced by phage display derived recombinant antibodies.

Applications

Immunofluorescence